



## Skin and Wound Product Information Sheet

<b>Biatain Ibu</b>	
<b>Classification</b>	<b>Analgesic Dressing: Ibuprofen</b>
<b>Key Points</b>	<ul style="list-style-type: none"> <li>• Polyurethane foam dressing which contains Ibuprofen (Ibu) 0.5gm/cm<sup>2</sup> homogeneously dispersed throughout the dressing.</li> <li>• Ibuprofen is released to the wound bed when in contact with wound exudate.</li> <li>• Dressing is designed to fully conform to a wound depth up to 2 cm due to the dressing's 3DFit Technology. This conformity allows the dressing to act as a primary cover dressing for a wound with a visible wound bed as it ensures no dead spaces between the wound bed and the dressing which could allow for the pooling of exudate.</li> </ul>
<b>Indications</b>	<ul style="list-style-type: none"> <li>• For use as a primary cover dressing for painful wounds.</li> </ul>
<b>Precautions</b>	<ul style="list-style-type: none"> <li>• Consult with Wound Clinician/NP/Physician prior to use.</li> <li>• Ensure that evaporating solutions (e.g., alcohol containing products) are completely dried off before dressing application.</li> </ul>
<b>Contraindications</b>	<ul style="list-style-type: none"> <li>• Do not use for clients with known sensitivity to polyurethane, Ibuprofen, acetylsalicylic acid or other NSAIDs.</li> <li>• Do not use on wound bed that is covered with more that 30% necrotic tissue.</li> <li>• Do not use with hydrogen peroxide or Dakin's Solution.</li> <li>• Do not use during treatment with radiation (x-ray, ultrasound, cancer radiation therapy).</li> <li>• Do not use on untreated clinically-infected wounds.</li> <li>• Do not use on children under 12 years of age except on the advice of a physician.</li> </ul>
<b>Formats &amp; Sizes</b>	<ul style="list-style-type: none"> <li>• Non-Adhesive Dressing                             <ul style="list-style-type: none"> <li>▪ 10 x 10 cm</li> <li>▪ 10 x 20 cm</li> <li>▪ 15 x 15 cm</li> <li>▪ 20 x 20 cm</li> </ul> </li> <li>• Soft hold Dressing                             <ul style="list-style-type: none"> <li>▪ 10 x 10 cm</li> <li>▪ 10 x 20 cm</li> </ul> </li> </ul> <div style="text-align: right; margin-top: 10px;">  </div>
Application Directions	Rationale
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
If required, apply a skin barrier to peri-wound skin as needed.	To protect the peri-wound skin from maceration and to improve the adhesion of the dressing or tape.
Choose appropriate dressing size; the dressing should extend at least 2.5cm beyond wound margin for effective absorption. Dressing may be cut if needed.	Incorrect sizing will adversely affect dressing function.
To Apply	
If there is minimal exudate, moisten the Biatain Ibu with a small amount of normal saline or apply a small amount of amorphous gel to the wound bed.	Ibuprofen is released to the wound bed when in contact with wound exudate.
Place the dressing directly over the wound; ensure the plain non-printed side is applied to the wound.	The outer side prevents strike through of exudate.
Secure with tape or wrap with kling/tape.	
To Remove	
Carefully lift the dressing off the wound.	To avoid trauma to the peri-wound skin.
Frequency of Dressing Change	
Will depend upon the amount of exudate. Dressing should be changed when the exudate spreads within 2.5cm of the edge of the dressing.	The absorbed drainage is clearly visible through the backing of the dressing.



## Skin and Wound Product Information Sheet

Application Directions	Rationale
Can be left on up to 7 days if pain is effectively managed; normal wear time is 2-3 day.	
Expected Outcome	
Wound pain is managed within 2-3 days.	
<b>For further information, please contact your Wound Clinician.</b>	